




Surgical Instruments **CATALOGUE** *(Part 3)*


*"2S Leverantörer Tand och
Sjukhus Instrument Handelsbolag"*






Ophthalmic Instrument

Barraquer Wire Speculum




IM-1001 (43101-03)
Small-Blade Length:10mm

Barraquer Heavy Wire Speculum




IM-1002 (43101-04)
Adult Size-Blade Length:14mm

Barraquer Speculum with Solid Blade




IM-1003
Blade Length:11mm




IM-1004
Blade Length:15mm

Kratz-Barraquer Heavy Wire Speculum with Open Blade

Easy Access to the Cataract Incision
with Phaco or I/A Tip




IM-1005
Small




IM-1006
Large

Wire Speculum




IM-1007
Temporal
Spring toward Nasal Side




IM-1008
Kratz-Barraquer
Blade Length:18mm

McIntyre Wire Speculum

Gently curved springs
V-shaped blades




IM-1009
Closed Wire:13mm Blades



IM-1010
Open Wire:13mm Blades

Nevyas Wire Speculum



IM-1011
Infant-Blade Length:4mm

Frankel Infant Lid Speculum

7mm

5mm

IM-1012
One Blade 5mm, One Blade 7mm


Lieberman Style Wire Speculum

Adjustable Mechanism




IM-1013,IM-1013T

Blades Length:15mm, K-wire, Standard




IM-1014,IM-1014T

Blades Length:15mm
with Traction Suture Anchor Posts




IM-1015,IM-1015T

Blades Length:10mm, K-wire, Small




IM-1016,IM-1016T

Blades Length:15mm, K-wire, Flattened




IM-1017,IM-1017T

Blades Length:11mm, Rounded



IM-1018,IM-1018T

Blades Length:13mm, Triple Post




IM-1019,IM-1019T

Blades Length:15mm, V-shaped

Lieberman Temporal Approach Speculum


Adjustable Mechanism on Nasal Side
Open Wire Blades



IM-1020,IM-1020T

Blades Length:15mm

Double-X Speculum



IM-1021
Blade 15mm, K-wire, Small



IM-1022
Blade 10mm, K-wire



IM-1023
Blade 15mm, V-wire




IM-1024
Blade 15mm, Solid Blades

Castroviejo Eye Speculum



IM-1025 (43132-08)
Blade 14mm, Adjustable Mechanism


Schott Eye Speculum



IM-1026
16mm wire blades

Murdoch Eye Speculum

With Locking Screw



IM-1027(43115-05)
Blades Length:12mm

IM-1028
Blades Length:15mm

CookEyeSpeculum

With LockingScrew




IM-1029
Blades Length:8mm

IM-1030
Blades Length:11mm

IM-1031
Blades Length:15mm

Lancaster Eye Speculum

Adjustable Mechanism
with Locking Nut



IM-1032(43129-07)
Blades Length:16mm

Lieberman Eye Speculum



IM-1033,IM-1033T
Solid Blades, Adjustable

Kraff Eye Speculum



IM-1034, IM-1034T
Nasal Solid Blades, Adjustable

Slade LASIK Speculum



IM-1035,IM-1035T
Blade Length:17mm

Kraff Nasal LASIK Speculum



IM-1036, IM-1036T
Blade Length:19mm

Slade Aspirating LASIK Speculum

Used with Suction Set
on Lower Syringe



IM-1037
Blade Length:17mm

Machat-Buratto Aspirating LASIK Speculum


Used with Suction Set
on Lower Syringe



IM-1038
Blade Length:15mm


Lieberman Aspirating Speculum

Used with Suction Set
on Lower Syringe



IM-1039
Blade Length:15mm

Flieringa Scleral Fixation Ring



IM-1061
Set of eight:15mm-22mm

Schepens Scleral Depressor

7mm Wide Tip




IM-1062
Small Thimble

IM-1063
Medium Thimble

IM-1064
Large Thimble

SchocketScleralDepressor

Doubleendedwith
Pocketclip



IM-1065

Jaeger Lid Plate



IM-1066(43150-11)

03

Stevenson Lacrimal Sac Retractor

Self-Retaining
Adjustable



IM-1081(43139-07)
3x3 Blunt Prongs

Stevenson Lacrimal Sac Retractor

Self-Retaining
Adjustable




IM-1082
Solid Blades with Serrated Edge

Jaffe Lid Retractor



IM-1083
Child, 9mm Blades(Pair)



IM-1084
Adult, 14mm Blades (Pair)

Agricola Lacrimal Sac Retractor



IM-1085
3x3 Sharp Prongs

Knapp Lacrimal Retractor



IM-1086
4-Prong, Blunt

IM-1087
4-Prong, Sharp

Desmarres Lid Retractor



IM-1088
Size-0,8mm

IM-1089
Size-1,10mm

IM-1090
Size-2,12mm

IM-1091
Size-3,14mm

Schepens Orbital Retractor



IM-1092

Walton Conjunctiva Retractor



IM-1093,IM-1093T
Blade Length:6mm wide,SolidThin

Helveston Tissue Retractor

Thin Curved Blade




IM-1094,IM-1094T
Blade Length:7mm wide

IM-1095,IM-1095T
Blade Length:9mm wide

IM-1096,IM-1096T
Blade Length:11mm wide,Concave

Castroviejo Caliper



IM-1097,IM-1097T(43583-07)
Straight, 20mm Spread

IM-1098,IM-1098T
Curved, 20mm Spread

Braunstein Caliper



IM-1099,IM-1099T
3.0mm & 3.5mm

IM-1100,IM-1100T
3.5mm & 4mm

Wright Fascia Needle



IM-1101
140mm Long, 1mm x 6mm Oval Hole

Bowman Malleable Lacrimal Probe

IM-1121
Size 0000 & 000


IM-1122
Size 00 & 0

IM-1123
Size 1 & 2

IM-1124
Size 3 & 4

IM-1125
Size 5 & 6

IM-1126
Size 7 & 8



IM-1127
Set of 6 Sizes

Infant Lacrimal Dilator



IM-1128, IM-1128T
Very Fine Blunt Tip

Muldoon Lacrimal Dilator



IM-1129, IM-1129T
Blunt Tip

Wilder Lacrimal Dilator



IM-1130, IM-1130T
Size-1, Fine Taper

IM-1131, IM-1131T
Size-2, Medium Taper

IM-1132, IM-1132T
Size-3, Heavy Taper

Castroviejo Lacrimal Dilator



IM-1133
Double-end fine and medium taper

Pigtail Probe



IM-1134
With Suture Holes

Harms Trabeculotomy Probe



IM-1135
With Suture Holes
Left

IM-1136
With Suture Holes
Right

Crawford Retrieving Hook



IM-1137

Beyer Rongeur



IM-1138
Single Action

Nasal Speculum




IM-1139
Infant

IM-1140
Child


IM-1141
Adult

Joseph Periosteal Elevator




IM-1142
Gently Curved Tip

Mallet Lead Filled




IM-1143

Troutman BladeBreaker & Holder



IM-1161, IM-1161T
Small

Castroviejo Blade Breaker & Holder



IM-1162, IM-1162T
Large

Scalpel Handles



IM-1163
Suitable for blades

Scalpel Handles



IM-1164
16.5cm

IM-1165
12cm

Scalpel Handles



IM-1166
10cm

IM-1167
13cm

Freer Periosteal Elevator



IM-1181
One End Blunt & Other End Semi Sharp

Stallard Dissector



IM-1182
Blunt

Rollet Rugine



IM-1183

West Lacrimal Chisel



IM-1184
Straight

Feer Lacrimal Chisel



IM-1185
Curved

Tooke Corneal Knife



IM-1186
Straight with Curved Cutting Edge

Corneal Dissector



IM-1187
Straight Blade, Angled

IM-1188
Curved Blade, Angled

Gill Corneal Knife

Curved Cutting Edge



IM-1189
Straight

IM-1190
Light Curve

DALK Lamellar dissector




IM-1191
Semi-sharp dissector helps remove the stromal fibers

Fukasaku "Hockey" Knife



IM-1192
Sharp 3mm Front Edge, Semi-sharp 6mm Bottom Edge

Foreign Body Spud




IM-1193 Golf Club
Angled, Lightly Round Tip

IM-1194 Francis
Angled, Flag Shaped Tip


IM-1195 Ellis
Lightly Curved, Round Tip

Spud & Gouge



IM-1196
Reversible Handle

Kelly Membrane Punch



IM-1201, IM-1201T
0.8mm Dia. Tip

Gass Corneoscleral Punch

Rotatable for
Convenient Positioning



IM-1202
1.5mm Dia. Tip


Fukasaku Micro Trabeculectomy Punch

Bullet Shaped Tip
0.3mm x 0.6mm Deep Bite



IM-1203T
0.6mm Dia. Head

Kerrison Rongeur




2/3

IM-1204
Size-1, 3.0mm Wide

IM-1205
Size-2, 4.0mm Wide

Citelli Rongeur



2/3

IM-1206
Size-2, 3.0mm Wide

Mc Pherson Iris Spatula



IM-1221
Malleable, 1mm Wide

Culler Iris Spatula

Malleable Blades



IM-1222
Malleable, 1mm Wide

IM-1223
Malleable ,2mm Wide

Wheeler double Ended Spatula



IM-1224
Lightly Curved

Elschnig Cyclodialysis Spatula



IM-1225
Tip Width:1mm

Barraquer Iris Spatula



IM-1226
Tip Width 0.50mm
Round, Smooth

IM-1227
Tip Width 0.25mm
Round, Smooth

Castroviejo Synechiaie Spatula




IM-1241
Double Ended, Semi-Sharp Tips, Tip Width: 0.50mm

Lewis & SNELLEN Lens Loop



IM-1243
Small



IM-1244
Large



IM-1245(43443-13)
Serrated

Kirby Lens Loop



IM-1246
Lens Loop and Muscle Hook

Gass Retinal Detachment Hook



IM-1247
With oval hole

Stevens Tenotomy Hook



IM-1248
Small,Blunt

Graefe Muscle Hook




IM-1249
Size-1

IM-1250
Size-2

IM-1251
Size-3

Jameson Muscle Hook



IM-1252
Small

IM-1253(43353-13)
Large

Green Muscle Hook



IM-1254(43355-13)


Fixation Hook



IM-1255
Double,Sharp,Small

IM-1256
Double,Sharp,Large

Meyerhoefer Chalazion Curette




IM-1266 Size-0 1.50mm Dia.	IM-1267 Size-1 1.75mm Dia.	IM-1268 Size-2 2.0mm Dia.
IM-1269 Size-3 2.50mm Dia.	IM-1270 Size-4 3.50mm Dia.	

Culler Lens Spoon



IM-1271
SterlingSilver

Bunge Evisceration Spoon



IM-1272 Small	IM-1273 Large
------------------	------------------

Wells Enucleation Spoon



IM-1274

Bonn Microlris Hook



IM-1281,IM-1281T

Tip Width:0.55mm,Very Delicate,Straight Blunt Tip

Hirshman Iris Hook



IM-1282,IM-1282T

Blunt,Angled

Sinsky II IOL Manipulating Hook



IM-1283,IM-1283T(43388-02)

Straight with Round Handle
0.2mm Dia. Blunt Tip

IM-1284,IM-1284T(43389-02)

Angled with Round Handle
0.2mm Dia. Blunt Tip

Jaffe Lens Manipulating Hook



IM-1285,IM-1285T

Straight with Flat Handle
0.15mm Dia. Blunt Tip

IM-1286,IM-1286T

Angled with Flat Handle
0.15mm Dia. Blunt Tip


Lewicky IOL Manipulating Hook



IM-1287,IM-1287T

Angled, 0.15mm Dia. Blunt Tip

Rentsch IOL Manipulating & Iris Hook



IM-1288, IM-1288T
Angled

Kuglen Iris Hook & Lens Manipulator



IM-1289, IM-1289T
Straight, Push-Pull Model

IM-1290, IM-1290T
Angled, Push-Pull Model

Kuglen Iris Hook & Manipulator



IM-1291, IM-1291T (43390-05)
Straight

Fenzl Lens Manipulating Hook



IM-1292, IM-1292T
Straight, 0.15mm Dia. T-Shaped Tip

IM-1293, IM-1293T
Angled, 0.15mm Dia. T-Shaped Tip

Maltzman-Fenzl Lens Manipulating Hook



IM-1294, IM-1294T
Angled, 0.15mm Dia. V-Shaped Tip

Jaffe-Maltzman Lens Manipulating Hook



IM-1295, IM-1295T
Angled, 0.15mm Dia. V-Shaped Tip, Flat Handle

Bechert Nucleus Rotator



IM-1296, IM-1296T
Angled 7mm, Blunt tip

IM-1297, IM-1297T
Angled 9mm, Blunt tip

Clayman Nucleus Rotator



IM-1298, IM-1298T
0.5mm wide flat spatula
with v-notched tip

Lester IOL Manipulator



IM-1299, IM-1299T
0.25mm dia. tip, straight

IM-1300, IM-1300T
0.25mm dia. tip, angled

Jaffe IOL Spatula



IM-1301, IM-1301T
Angle to Tip: 10mm, Tip Width: 0.50mm, Flat Smooth Tip

Escaf Nucleus Manipulator



IM-1302,IM-1302T
50° Angled, Flattened Tip

Knolle Nucleus Spatula



IM-1303,IM-1303T
0.5mm Wide, Flat

Koch Nucleus Spatula



IM-1304,IM-1304T
Gently Curved, Smooth

IM-1305,IM-1305T
Gently Curved, Notches



Escaf Nucleus Sustainer



IM-1306,IM-1306T
Smooth Ball Tip

Graether Collar Button Micro Iris Retractor



IM-1307,IM-1307T
Straight



IM-1308,IM-1308T
Angled

Lieberman Micro Finger & Spatula



IM-1309,IM-1309T

Lieberman Micro Finger and Spatula





IM-1310,IM-1310T
For Use in Right Hand

IM-1311,IM-1311T
For Use in Left Hand

Fukasaku-LiebermanPhacoSpatula

Small pupil snapper hook
and the Micro Finger



IM-1312,IM-1312T
For Use in Right Hand

IM-1313,IM-1313T
For Use in Left Hand

Grayson Nucleus Manipulator



IM-1314,IM-1314T
1.1mm Long Tip, Flattened

Drysdale Nucleus Manipulator



IM-1315,IM-1315T
Rounded Edges

IM-1316,IM-1316T
Defined Interior Edge

Tennant Nucleus Manipulator



IM-1341,IM-1341T
Tip Dia. 0.6mm, To Split and Rotate Nucleus

Nagahara Phaco Chopper



IM-1342,IM-1342T
Cutting Edge 90°, Cutting Length:1.5mm

Nagahara Phaco Chopper



IM-1343,IM-1343T
Cutting Edge 60°, Cutting Length:1.5mm

Koch-Nagahara Spatula & Chopper



IM-1344,IM-1344T
Cutting Edge 60°, Cutting Length:1.5mm

Kammann Prechopper



IM-1345,IM-1345T
Angled, Cutting Edge 2.2mm

IM-1346,IM-1346T
Straight, Cutting Edge 2.2mm

Friedman-Koch Universal Chopper & Spatula



IM-1347,IM-1347T
Double Cutting Edges
Cutting Edge(Chopper):1.5mm, Tip Width (Manipulator):0.8mm

OlsonPhacoChopper



IM-1348,IM-1348T
Cutting Edge60,°Cutting Length:2mm

Rosen Phaco Splitter



IM-1349,IM-1349T
WedgeShaped

Nordan Capsule Elevator



IM-1350,IM-1350T
Tip Dia. 1.5mm

Nagahara Bi-Manual Irrigating Chopper




IM-1351T
Round Tip, 20Gauge Tube

Nagahara Bi-Manual Irrigating Carate Chopper



IM-1352T
Sharp Tip, 20 Gauge Tube

C MICS Phaco Pre- hopper



IM-1353,IM-1353T
Two Defined Edges

IM-1354,IM-1354T
Curved with Blunt Tip
and Sharp Inferior Edge

Phaco Chopper



IM-1355,IM-1355T
Ball Tip, 0.7mm Inferior Wedge


Nagahara Phaco Chopper



IM-1356,IM-1356T
Rounded, Blunt Tip, 0.9mm Inferior Wedge

Chang-Seibel Phaco Chopper

Acombination of chopper with the
Chang Micro Finger and end
the Seibel Vertical Safety Chopper



IM-1357,IM-1357T
For Use in Right Hand

IM-1358,IM-1358T
For Use in Left Hand

Milder Phaco Chopper


60° off set in ferior wedge
for horizonta IMhopping and a
beveled distal tip for vertical chopping



IM-1359,IM-1359T
For Use in Left Hand

MilderPhacoChopper


60° off set in ferior wedge
for horizontal chopping and a
beveled distal tip for vertical chopping
12mm long spatula with disc-shaped tip



IM-1360,IM-1360T
For Use in Right Hand

IM-1361,IM-1361T
For Use in Left Hand

Masket Phaco Spatula



IM-1362,IM-1362T
12mm Long Round Spatula, Disk-shaped Tip

Chang Phaco Chopper

0.9mm Long Quick Chop Tip



IM-1363,IM-1363T
For Use in Right Hand



IM-1364,IM-1364T
For Use in Left Hand

Minardi Phaco Chopper



IM-1365,IM-1365T
0.7mm Long Quick Chop Tip

Hwang Phaco Chopper



IM-1366,IM-1366T
1mm Long Vertical Chopper and a Chisel Tipped Spatula

Agarwal Karate Phaco Chopper





IM-1367,IM-1367T
Semi Sharp, Pointed Tip

Chang Micro Finger Chopper



IM-1368,IM-1368T
Curved Finger Shaped Tip

Warren Canine Phaco Chopper




IM-1369,IM-1369T
2mm Pointed Tip, 3.5mm Long Micro Finger

Brierley Nucleus Splitter



IM-1370,IM-1370T
Incision Sealing Shaft

Dodick-Kammann Nucleus Chopper



IM-1371,IM-1371T	IM-1372,IM-1372T	IM-1373,IM-1373T	IM-1374,IM-1374T
1.5mm Long Tip, Straight	1.5mm Long Tip, Angled	2mm Long Tip, Straight	2mm Long Tip, Angled

Koch Nucleus Chopper



IM-1375,IM-1375T
1.6mm Long Blunted Tip

Escaf Nucleus Divider




IM-1376,IM-1376T
1.25mm Long Wedge Shaped Tip

Neuhann Nucleus Divider



IM-1377,IM-1377T
Strongly Curved Spatulated Tip

DALK Corneal Dissector



IM-1391,IM-1391T
Flat,12mm Long, Beveled

IM-1392,IM-1392T
Triangular Shaft

IM-1393,IM-1393T
1mm Wide, 8mm Long, Flattened

Gimbel Nucleus Spatula



IM-1394,IM-1394T
0.25mm Dia. ,Smooth Angled 30°

Girard Synechiaie Spatula



IM-1395,IM-1395T
0.5mm Wide Tip, Semi Sharp

Obstbaum Synechiaie Spatula




IM-1396,IM-1396T
0.25mm Wide Tip, Semi Sharp

Benedet to Intra Lase Spatula




IM-1397,IM-1397T
10mm Long, Flattened with Beveled Notches

Moran LASIK Enhancement Spatula



IM-1398,IM-1398T
A Modified Paton Spatula

Hersh Intra Lase Flap Spatula



IM-1399,IM-1399T
6mm & 10mm, Beveled Tips

Binder Intra Lase Flap Elevator



IM-1400,IM-1400T
1mm Wide, 10mm Long, Semi-sharp Edges

Fukasaku LASIK Spatula




IM-1401,IM-1401T
Rounded Interior, Beveled Edges

Guimaraes LASIK Spatula



IM-1402,IM-1402T
Gently Curved, Bi-convex, Semi-sharp Edges

Katzen LASIK Flap “Unzipper”



IM-1403,IM-1403T
Bi-convex Spatula with Beveled Notches

Maddox LASIK Spatula



IM-1404,IM-1404T
Flat & Cylindrical Spatula


Rowen Combination LASIK Spatula



IM-1405,IM-1405T

Endothelial Stripper


1mm Wide, Thin Beveled Tip



IM-1406,IM-1406T
Smooth Surface

IM-1407,IM-1407T
Textured Surface


ThorlaksonZ-LASIK Flap Lifter



IM-1408,IM-1408T
Beveled Notches, Long Handle

IM-1409,IM-1409T
Beveled Notches, Short Handle

Diamond Knives



IM-1410T
Straight/Single Edge
FH-MM6, 0x1, 0x0, 17mm/30°



IM-1411T
Tri-facet
FH-MDS5, 1x1, 0x0, 20mm/45°/30°



IM-1412T
Straight/Single Edge
FH-MM4, 0x0, 7x0, 17mm/45°



IM-1413T
Tunnel
FH-RS19-facets 6, 0x2, 0x0, 17mmabg



IM-1414T
ClearCorneaMICS
FH-CC6, 0x2, 2x0, 20mm/90° abg



IM-1415T
Clear Cornea 90°
FH-CC6, 0x3, 0x0, 20mm/90° abg




IM-1416T
Tapered Spear MICS
FH-CCtr6, 0x2, 4/2, 70x0, 20mm/90° abg



IM-1417T
Clear Cornea 60°
FH-ZALDIVAR6, 5x0, 55/1, 0X0, 17mmabg

3D Diamond Knives



IM-1418T
FH-CC3D6, 0x2, 20x0, 17mm/80° abg



IM-1419T
FH-CC3D6, 0x2, 40x0, 20mm/80° abg




IM-1420T
FH-CC3D6, 0x2, 50x0, 20mm/80° abg



IM-1421T
FH-CC3D6, 0x2, 65x0, 17mm/80° abg

Diamond Knife with Glaucoma Surgery




IM-1422T

Step Knife

Steps 0.15/0.2/0.25/0.3/0.37/0.5/4.0

ST-3 Schnitt6, 0x1, 0x0, 20mm/25°

Diamond Knife with LRI Surgery




IM-1423T

With 3 Settings: 0.5, 0.55, 0.6

ST-3 Schnitt5, 1x0, 4/1, 0x0, 10mm/20°

Step Diamond Knife



IM-1424T

Tri-facet

Step sat 0/0.3/0.4/0.5/0.6/0.8/1.0/2.5mm

ST-3 MDS5, 6x1, 0x0, 17mm

Hofmann-Thornton Globe Fixation Ring



IM-1431,IM-1431T
13mm Dia. Multiple Blunt Teeth

IM-1432,IM-1432T
16mm Dia. Multiple Blunt Teeth

Fine-Thornton Swivel Fixation Ring



IM-1433,IM-1433T
13mm Dia. Blunt Teeth

IM-1434,IM-1434T
16mm Dia. Blunt Teeth

Gimbel Stabilizing Ring



IM-1435,IM-1435T
13mm Dia. ,Blunt Teethon Both Sides

Thornton Globe Fixation Ring



IM-1436,IM-1436T
16mm Dia. Multiple Blunt Teeth

Lu-Mendez LRI Guide and Fixation Rings



IM-1437,IM-1437T
Delicate Blunt Teeth, Marks 10° Increments, 11mm I.D.

Mendez Degree Gauge



IM-1438,IM-1438T
12mm I.D. Marked in 10° Increments

Bores Axis Marker



IM-1439,IM-1439T

Camellin Style LASEK AIMohol Well



IM-1440,IM-1440T
8.5mm Dia.

IM-1441,IM-1441T
9.5mmDia.

McDonald Optic Zone Marker



IM-1442,IM-1442T
6.5mm & 7mm Dia., Double End, CrossHairs

Anis Radial Marker



IM-1443,IM-1443T
4mm Dia. 8 Blades

IM-1444,IM-1444T
4mm Dia. 12 Blades

IM-1445,IM-1445T
4mm Dia.16 Blades

Toric Marker



IM-1446,IM-1446T
Half-ring for mark of meredian 12mm

Intra-Op Toric Axis Marker



IM-14 47T
Markings every 5 Degrees

Kellen Capsulorhexis Marker



IM -1448, IM-14 48T
6mm Dia.

IM -1449, IM-1449T
5.5mm Dia.

Osher Neuman Corneal Marker



IM -1450, IM-1 450T
4 Blades

IM -1451, IM-1 451T
6 Blades

IM -1452, IM-1 452T
8 Blades

IM -1453, IM-1
16Blades

Hoffer Corneal Optical Zone Markers

FineCrossHairs



⊕ IM -1454, IM-1 454T
3.0mm

⊕ IM -1455, IM-1 455T
3.5mm

⊕ IM -1456, IM-1 456T
4.0mm

⊕ IM -1457,
4.5mm

⊕ IM -1458, IM-1 458T
5.0mm

⊕ IM -1459, IM-1 459T
5.5mm

⊕ IM -1460, IM-1 460T
6.0mm

⊕ IM -1461,
6.5mm

⊕ IM -1462, IM-1 462T
7.0mm

⊕ IM -1463, IM-1 463T
7.5mm

⊕ IM -1464, IM-1 464T
8.0mm

Maloney Intra-Operative Keratometer



IM -1465, IM-1465T
Measurement of Astigmatism

Iris Scissors



IM-1469(09340-10)

Straight, Delicate, Pointed Tips 10cm

IM-1470(09341-10)

Curved, Delicate, Pointed Tips 10cm

Iris Scissors




IM-1471(09340-11)

Straight, Standard, Pointed Tips 11cm

IM-1472(09341-11)

Curved, Standard, Pointed Tips 11cm


Stevens Tenotomy Scissors



IM-1473

Needle Points, Lightly Curved

Stitch Scissors



IM-1474 (09351-09)
Ribbon type, Blunt Tips, Straight

IM-1475
Ribbon type, Blunt Tips, Curved

Knapp Strabismus Scissors



IM-1476
Straight, Blunt Tips 11cm

IM-1477
Curved, Blunt Tips 11cm

Enucleation Scissors



IM-1478
Blunt tips, Medium Curve
14.5cm Long

IM-1479
Blunt tips, Strong Curve
14.5cm Long

Westcott Tenotomy Scissors



IM-1480
Blunt

IM-1481
Pointed

Noyes Scissors



IM-1482
Straight

Barraquer Iris Scissors



IM-1483
Blunt tips, 7mm on glades

Castroviejo Keratoplasty Scissors



IM-1484
Angle to Side: Blunt Tip

Westcott Type Scissors



IM-1485T
Straight, Blade Length 6mm
Flat Titanium Handle Double Action
Stainless Steel Tips

IM-1486T
Curved, Blade Length 6mm
Flat Titanium Handle Double Action
Stainless Steel Tips

Mc Pherson-Westcott Conjunctival Scissors



IM-1487
Blunt tips. Small blades, Curved

Jaffe Stitch Scissors



IM-1488
Very sharp pointed tips, Medium blades, Curved

Westcott Tenotomy Scissors



IM-1489
Blunt tips, Medium blades, Curved

Iris Scissors



IM-1490
Sharp tip, Small blades, Straight

IM-1491
Sharp tip, Small blades, Curved

Fine Stitch Scissors



IM-1492
Sharptips, Medium blades
Straight

IM-1493
Sharptips, Medium blades
Curved

Castroviejo Universal Corneal Scissors



IM-1494
Curved, Blunt Tip
Small blades

IM-1495
Curved, Blunt Tip
Medium blades

IM-1496
Curved, Blunt Tip
Large blades

IM-1497
Strongly Curved, Blunt Tip
Large blades

Troutman-Castroviejo Corneoscleral Scissors



IM-1498
With top, Small blades, Right

IM-1499
With stop, Small blades, Left

Katzin Corneal Transplant Scissors



IM-1500
Strong curve, Medium blades
Right

IM-1501
Strong curve, Medium blades
Left

Castroviejo Corneoscleral Scissors



IM-1502
Small blades, Right

IM-1503
Small blades, Left

Rich Corneal Transplant Scissors

Extremely Delicate
Strongly Curved Blades




IM-1504
Right

IM-1505
Left

Terry Micro Corneal Scissors

Round Handle



IM-1506
Right

IM-1507
Left

Tibolt Punctal Scissors



IM-1508
One conically shaped blade with grooves
and textured surface

DALK Micro Scissors Left+Right Set



IM-1509
The lower blade of the scissors is thinner
and 1mm longer than the upper blade

Vannas Capsulotomy Scissors



IM-1510



IM-1511



IM-1512



IM-1513



IM-1514

IM-1510,IM-1510T
Straight, Blade Length:7mm

IM-1511,IM-1511T
Curved, Blade Length:7mm

IM-1512,IM-1512T
Angled-Forward, Blade Length:7mm

IM-1513,IM-1513T
Angle to Side,Blade Length:7mm

McPhersonVannas
IM-1514,IM-1514T
Sharp Tips,Blade Length:8mm

Capsulotomy Scissors



IM-1515



IM-1516



IM-1517



IM-1518



IM-1519



IM-1520



IM-1521



IM-1522



IM-1523

7mm extra thin long blades, Sharp tips

IM-1515,IM-1515T
Vannas Scissors, Straight

IM-1516,IM-1516T
Vannas Scissors, Curved

10mm extra thin long blades, Sharp tips

IM-1518,IM-1518T
Vannas Scissors, Straight

IM-1520,IM-1520T
Stern-Gills,Angledforward

10mm long angled blades, Curved tips
IM-1522,IM-1522T
Gills, Sharp points

IM-1517,IM-1517T
Gills-Vannas, Angled forward

IM-1519,IM-1519T
Vannas Scissors, Curved

IM-1521,IM-1521T
Uribe-Stern, Curved, Angledforward

IM-1523,IM-1523T
Uthoff-Gills, Blunt points

Vitreo-Retinal Scissors



Angled 45°



IM-1524T
20 Gauge
IM-1525T
23 Gauge

Curved (R12mm)



IM-1526T
20 Gauge
IM-1527T
23 Gauge

Straight Blunt Tips



IM-1528T
20 Gauges
IM-1529T
23 Gauge

Straight Pointed Tips



IM-1530T
20 Gauges
IM-1531T
23 Gauge

Vertical 40°



IM-1532T
20 Gauge
IM-1533T
23 Gauge

Vertical 80°



IM-1534T
20 Gauges
IM-1535T
23 Gauge

Vertical, Lower Blade with Bead Tip



IM-1536T
20 Gauge
IM-1537T
23 Gauge

Alio MICS Scissors



IM-1538T
23 Gauge, Curved Shaft Length 28mm
with Horizontal Micro Blades, Tips Length 2.7mm

Dalk Air Cannula



IM-1560
27 Gauge, Flattened

Rycroft Air Injection Cannula



IM-1561
30 Gauge, Angle to Tip:5mm

IM-1562
27 Gauge, Angle to Tip:5mm

Langs Lacrimal Cannula



IM-1563
Curved, Long

Langs Lacrimal Cannula



IM-1564
Curved, Short


Harrison Butler Lacrimal Cannula



IM-1565

Binkhorst Irrigation Cannula

Smooth Blunt Tips
22 Gauge




IM-1566
Straight

IM-1567
Left

IM-1568
Right

McIntyre Anterior Chamber Cannula



IM-1569 25 Gauge, Angled	IM-1570 23 Gauge, Straight	IM-1571 23 Gauge, Angled	IM-1572 26 Gauge, Straight
IM-1573 26 Gauge, Angled	IM-1574 30 Gauge, Straight	IM-1575 30 Gauge, Angled	


Bishop-Harmon Irrigation Cannula



IM-1576
Angle to Tips: 7mm, 19 Gauge

Simcoe Irrigation-Aspirating Cannula

for 120° Clock Aspiration, 23 Gauge
0.3mm Aspirating Side Port
with Silicone Tubing and Adapter



IM-1577 Right	IM-1578 Left
------------------	-----------------

Jensen Posterior Capsule Polisher



IM-1579 Sand Blasted Tip, 23 Gauge	IM-1580 Sand Blasted Tip, 25 Gauge
---------------------------------------	---------------------------------------

Capsule Polishing Curette


Semi-Sharp Edges for Anterior and Posterior
Capsule Polishing Irrigation in to Ring



IM-1581
25 Gauge

Gills Aspirating Cannula

Smooth blunt tips,
25 Gauge, Angled



IM-1582 5mm Tip IM-1583 7mm Tip IM-1584 8mm Tip IM-1585 10mm Tip

Simcoe Cortex Extractor

0.3mm Port Size, for Side Port
Cortex Removal



IM-1586 13mm, 23 Gauge IM-1587 21mm, 23 Gauge

Knolle-Pearce Irrigating Vectis




IM-1588 23 Gauge

Anis Irrigating Vectis



IM-1589 23 Gauge

Irrigating-Aspirating Cannula




IM-1590 Simcoe IM-1591 Reverse Simcoe

De La O Aspirating Cannula

IM-1592
0.3mm side port, Angled at 12mm, 25 Gauge

Morrison Aspirating Cannula

IM-1593
Smooth blunt tip, Angled at 7mm, 23 Gauge

Anis Cortex Aspirating Cannula

IM-1594
0.4mm side port, 23 Gauge, Straight

IM-1595
0.4mm side port, 23 Gauge, Angled

Wilbrandt-Binkhorst Cannula

IM-1596
Smooth blunt tip,
23 Gauge, Straight

IM-1597
Smooth blunt tip,
23 Gauge, Angled left

IM-1598
Smooth blunt tip,
23 Gauge, Angled right

McIntyre-Binkhorst Cannula

IM-1599
Smooth blunt tip,
26 Gauge, Straight

IM-1600
Smooth blunt tip,
26 Gauge, Angled left

IM-1601
Smooth blunt tip,
26 Gauge, Angled right

Bishop-Harmon Anterior Chamber Irrigator

Includes 19 gauge angled cannula, adapter
and blue silicon bulb

IM-1602
Complete

IM-1603
Siliconbulb

IM-1604
Adapter

IM-1605
19 gauge cannula

IM-1606
21 gauge cannula

Nucleus Hydrodis section Cannula



IM-1607
27 Gauge, Flat Angled Tip 8mm

IM-1608
25 Gauge, Flat Angled Tip 8mm

Kellan Hydrodis section Cannula



IM-1609
25 Gauge, Curved, Vertical Flattened Tip

Kellan Hydrodis section Cannula



IM-1610
25 Gauge, Curved Bevel Tip

LASIK Double Arm Cannula



IM-1611
25 Gauge, 8 Ports

IM-1612
Curved, 25 Gauge, 8 Ports

Guell LASIK Irrigation Cannula



IM-1613
25 Gauge, 7 Ports

Buratto LASIK Irrigation Cannula



IM-1614
25 Gauge, 3 Ports

Feaster Hydrodis section Cannula



IM-1615
25 Gauge, Wedge Shaped Tip

Chang Hydrodis section Cannula

27 Gauge, Flat 90° Angled Flat




IM-1616
Standard Tip

IM-1617
Long Tip

Simcoe Irrigating Lens Nucleus Loop

25 Gauge, Finely Serrated,
Strong Curve




IM-1618
Right

IM-1619
Left

McIntyre Lacrimal Cannula

23 Gauge, Dual 0.3mm Side Port





IM-1620
Straight

IM-1621
Curved

Lacrimal Cannula


23 Gauge, Malleable Tip



IM-1622
Straight

IM-1623
Curved


Heslin Infusion Cannula



IM-1624
25 Gauge


IM-1625
23 Gauge

Infusion Cannula



IM-1626
22 Gauge, Self Retaining Tip

Brierley Chamber Maintainer




IM-1627
20 Gauge, Tapered Oval Tip with Self Retaining Fins

Lewicky Chamber Maintainer



IM-1628
20 Gauge, Threaded Tip

Dahan-Blumenthal Anterior Chamber Maintainer



IM-1629
20 Gauge, Tapered Oval Tip

Buratto Bi-Manual Irrigation Hand piece

Open End



IM-1641T
21 Gauge, Straight

IM-1642T
23 Gauge, Curved

Buratto Bi-Manual Irrigation Hand piece

DualPort



IM-1643T
21 Gauge, Straight
Port Size 0.5mm

IM-1644T
21 Gauge, Curved
Port Size 0.5mm

IM-1645T
23 Gauge, Straight
Port Size 0.4mm

IM-1646T
23 Gauge, Curved
Port Size 0.4mm

Buratto Bi-Manual irrigating Handpiece



IM-1647T
21 Gauge, Curved
Port Size 0.4mm

IM-1648T
23 Gauge, Curved
Port Size 0.3mm

Buratto Bi-Manual Aspiration Handpiece



IM-1649T
21 Gauge, Curved
Port Size 0.4mm

IM-1650T
23 Gauge, Curved
Port Size 0.3mm

Buratto Bi-Manual Aspiration Handpiece



IM-1651T
21 Gauge, Port Size 0.35mm
Curved with Capsule Polisher

Gills-Colibri Forceps



IM-1661,IM-1661T
Very Fine Pointed Tips

Colibri Forceps



IM-1662,IM-1662T(43885-07)
Very Delicate, 1x2 Teeth, 0.12mm

Harms-Colibri Forceps

With Tying Platform



IM-1663,IM-1663T
Very Delicate, 1x2 Teeth, 0.12mm

Pierse-Colibri Corneal Forceps



IM-1664,IM-1664T
Very Delicate

Castroviejo-Colibri Forceps

Very Delicate, 1x2 Teeth, 0.12mm



IM-1665,IM-1665T
Without Tying Platform

IM-1666,IM-1666T
With Tying Platform


Polack Double Corneal Suturing Forceps



IM-1667,IM-1667T
1x2 Teeth, 0.1mm, Very Delicate

Bonn Forceps

With Tying Platform




IM -1668, IM-1668T
Short Handle, 1x2 Teeth, 0.12mm

IM -1669, IM-1669T
Long Handle, 1x2 Teeth, 0.12mm

McPherson Corneal Forceps

With Tying Platform



IM -1670, IM-1670T
Straight, 1x2 Teeth, 0.4mm

IM -1671, IM-1671T
Angled, 1x2 Teeth, 0.4mm

Kelman-McPherson Forceps



IM -1672, IM-1672T
Angled, 1x2 Teeth, 0.3mm, 8mm Long Tying Platform

Hoffer-McPherson Forceps



IM -1673, IM-1673T
Angled, Very Delicate, 1x2 Teeth, 0.12mm, 11 mm Long Tying Platform

Castroviejo Suturing Forceps


1x2 Teeth with Tying Platform



IM -1674, IM-1674T	IM -1675, IM-1675T	IM -1676, IM-1676T	IM -1677, IM-1677T	IM -1678,
0.12 mm	0.3mm	0.5mm	0.9mm	1.5mm

Castroviejo Suturing Forceps


Angled, 1x2 Teeth, With Tying Platform



IM -1679, IM-1679T
0.3mm

IM -1680, IM-1680T
0.5mm

Hoskin Corneal Forceps



IM-1681,IM-1681T
Straight, FineGrooved Tips

Hoskin Corneal Forceps



IM-1682,IM-1682T
Straight ,Micro Grooved Tips


Hoskin Corneal Forceps



IM-1683,IM-1683T
Straight, Colibri Style, Micro Grooved Tips

Mc Pherson Tying Forceps

5mm Long Smooth Tips
for 8-0 to 11-0 Sutures



IM-1684,IM-1684T
Straight

IM-1685,IM-1685T
Angled

Kelman-Mc Pherson Forceps



IM-1686,IM-1686T
Very Delicate, 7.5mm Long Smooth Tips

Bechert Forceps



IM-1687,IM-1687T
Angled, Very Delicate, 10mm Long Smooth Tips

Gaskin Fragment Forceps


11mm Long Smooth Tips
Very Delicate



IM-1688,IM-1688T
Angled

Tennant Tying Forceps

Extra Delicate, Smooth Tips
for 9-0 to 11-0 Sutures
6mm Long Platform




IM-1689,IM-1689T
Straight

IM-1690,IM-1690T
Curved

Harms Tying Forceps

for 7-0 to 9-0 Sutures
Smooth Tips



IM-1691,IM-1691T
Straight

IM-1692,IM-1692T
Curved

Jaffe Utility Forceps

for 7-0 to 9-0 Sutures



IM-1693,IM-1693T
Very delicate smooth jaws, Curved

Foreign Body Forceps




IM-1694,IM-1694T
Straight, Sharp Pointed Tips

Bonaccol to Utility & Splitter Forceps




IM-1695,IM-1695T
Straight, Very Delicate, Criss-Cross Serrated Tips

Mendez Multi-Purpose LASIK Forceps



IM-1696,IM-1696T
Smooth Grasping Surfaces, Spatulated Jaws and Tips

Utrata Capsulorhexis Forceps




IM-1697,IM-1697T
Very Delicate, Short Handle, Triangular Grasping Tips

Masket Capsulorhexis Forceps




IM-1698,IM-1698T
Very Delicate, Short Handle, Triangular Grasping Tips

Kraff Capsule Tag Forceps



IM-1699,IM-1699T
11mm Long Shanks, 1.5mm Long Tips

Giannetti MICS Capsulorhexis Forceps



IM-1700,IM-1700T
Ultra Delicate Shanks Tips, Maximum Jaw Opening 2mm

Utrata Capsularhexis Forceps



IM-1701T
Delicate Straight Shaft 11.5mm Long Shank
Round Handle 8mm Dia.

Utrata Capsularhexis Forceps



IM-1702T
Delicate Curved Shaft 11.5mm Long Shank
Round Handle 8mm Dia.

Utrata Capsularhexis Forceps



IM-1703T
Delicate Curved Shaft 10mm Long Shank
Round Handle 8mm Dia.

Utrata Capsularhexis Forceps



IM-1704
Extra Delicate Curved Shaft 11mm Long Shank
Round Handle 6mm Dia.

Dodick Nucleus Cracker



IM-1705,IM-1705T
Cross Action

Kraff Capsule Tag Forceps



IM-1706,IM-1706T
11mm Long Tips, 1.5mm Long Smooth Platform

Utrata Capsulorhexis Forceps



IM-1707,IM-1707T
Very Delicate Grasping Tips and Extremely Thin, 11mm Long Shanks

Utrata Capsulorhexis Forceps



IM-1708,IM-1708T
Round Handle, Very Delicate Grasping Tips and Extremely Thin 11mm Long Shanks

Masket Capsulorhexis Forceps



IM-1709,IM-1709T
Very Delicate Grasping Tips and Extremely Thin 11mm Long Shanks

Inamura Capsulorhexis Forceps

IncisionSealingDesign



IM-1710,IM-1710T
Short for Corneal Incisions

IM-1711,IM-1711T
Long for Scleral Incisions

MicroIncision Capsulorhexis Forcep


Squeeze Handle



IM-1712T
23 Gauge, 15mm Shaft Length

IM-1713T
23 Gauge, 25mm Shaft Length


Inamura Round Handle Capsularhexis Forceps for Scleral Tunnel Incision



IM-1714T

Curved Shaft, Tip to Joint 11mm
45° Angled Tip, Cross Action 1.5mm Wid that Box Joint


Inamura Round Handle Capsularhexis Forceps for Scleral Tunnel Incision



IM-1715T

Curved Shaft, Tip to Joint 10mm
45° Angled Tip, Cross Action 1.5mm Wid that Box Joint


Inamura Capsularhexis Forceps



IM-1716T

Spring Action, 45° Angled Tip
Curved Shaft, Tip to Joint 11mm

Inamura Capsularhexis Forceps



IM-1717T

Spring Action, 45° Angled Tip
Curved Shaft, Tip to Joint 8.5mm

Akahoshi Micro Prechopper Forceps



IM-1718T
Tips Height 1.3mm, Tips Open to 2.5mm
for Incision 1.7mm or more

Akahoshi Micro Prechopper Forceps



IM-1719T
Tips Height 1.5mm, Tips Open to 2.5mm
for Incision 1.7mm or more

Akahoshi Micro Prechopper Forceps



IM-1720T
Tips Height 1.6mm, Tips Open to 2.5mm
for Incision 1.7mm or more

Paddle Prechopper Forceps



IM-1721
Angled Shaft 45°, Tip to Angle Length 13.5mm
Tips Height 1.6mm, Tips Open to 2.5mm
for Incision 1.7mm or more

Fechtner Conjunctiva Forceps



IM-1722,IM-1722T

Very Delicate, Ring Shaped Tip with Tying Platform

Nugent Utility Forceps



IM-1723,IM-1723T

0.3mm Blunt Tips, Cross-action Mechanism

Stern-Conjunctiva Fixation Forceps




IM-1724,IM-1724T

1x2 Teeth 0.5mm With Lock

Moody Fixation Forceps

1x2 Teeth 0.5mm With Lock



IM-1725,IM-1725T

Left

IM-1726,IM-1726T

Right

Buratto LASIK Flap Forceps



IM-1727,IM-1727T

Disc Shaped,Serrated Jaws

Vitreo-Retinal Forcep

SqueezeHandle



Straight Forceps
with 1x2 Teeth



IM-1728T
20 Gauge

IM-1729T
23 Gauge

Serrated Crocodile Forceps
2mm & 4 mm Long Tips



IM-1730T
20 Gauge

IM-1731T
23 Gauge

Side Gripping Forceps
Smooth Round Jaws
2.4mm



IM-1732T
20 Gauge

IM-1733T
23 Gauge

Gripping Forceps,
2.5mm Smooth,
Straight



IM-1734T
20 Gauge

IM-1735T
23 Gauge

Membrane Peeler,
Smooth, 45°



IM-1736T
20 Gauge

IM-1737T
23 Gauge

Membrane Peeler,
Grooved, 45°



IM-1738T
20 Gauge

IM-1739T
23 Gauge

Straight and
End Gripping Jaws



IM-1740T
20 Gauge

IM-1741T
23 Gauge

Gripping Forceps,
2.0mm Smooth
Rounded Tips



IM-1742T
20 Gauge

IM-1743T
23 Gauge

ILM Forceps



IM-1744T
20 Gauge

IM-1745T
23 Gauge


Foreign Body Forceps,
Three Claws



IM-1746T
20 Gauge

IM-1747T
23 Gauge

Bishop-Harmon Tissue Forceps



IM-1748, IM-1748T
Delicate, 1x2 Teeth, 0.5mm

IM-1749, IM-1749T
Standard, 1x2 Teeth, 0.8mm


Bishop-Harmon Dressing Forceps



IM-1750, IM-1750T
Serrated Tips, 0.5mm


IM-1751, IM-1751T
Serrated Fine, 0.5mm

Sauer Fixation Forceps



IM-1752, IM-1752T
1x2 Teeth


Manhattan Fixation Forceps



IM-1753, IM-1753T
1x2 Angled Teeth with Tying Platform

Graefe Dressing Forceps

Serrated, Delicate




IM-1754
Straight

IM-1755
Curved

Graefe Tissue Forceps

1x2 Teeth ,Delicate



IM-1756
Straight

IM-1757
Curved

Elsching Scleral Forceps



IM-1758
Delicate

Remky Superior Rectus Forceps



IM-1759

Graefe Fixation Forceps



IM-1760
4.5mm Wide

Nugent Utility Forceps



IM-1761
Smooth Jaws

IM-1762
Serrated Jaws

Bonaccol to Conjunctiva Forceps





IM-1763
1.2mm wide, Serrated jaws

IM-1764
1.7mm wide, Serrated jaws

Adson

125mm Long, Straight



1x2 teeth, delicate

IM-1765
DressingForceps,
Serrated Jaws,
Delicate, 0.9mm

IM-1766
DressingForceps,
Serrated Jaws,
Standard, 1.3mm

IM-1767
TissueForceps,
1x2 Teeth,
Delicate, 0.9mm

IM-1768
TissueForceps,
1x2 Teeth,
Standard, 1.3mm

BayonetDressingForceps



IM-1769
146mm Long, Serrated Jaws

Green Fixation Forceps



IM -1770 (44076-05)
5mm Wide, Multiple Teeth

IM -1771
10mm Wide, Multiple Teeth

Barraquer Cilia Forceps




IM-1772
Lightly Angled, Smooth

Cilia Forceps



IM -1773
Slanted

Jewelers Forceps



IM -1774, IM-1774T
Standard, No.1

Jewelers Forceps



IM -1775, IM-1775T
Delicate, Short, No.3C

Jewelers Forceps



IM -1776, IM-1776T
Delicate, No.3

Jewelers Forceps



IM-1777,IM-1777T
Delicate, No.4

Jewelers Forceps



IM-1778,IM-1778T
Extra Delicate, No.5

Jewelers Forceps



IM-1779,IM-1779T
Extra Delicate, Angled, No.5

Jewelers Forceps



IM-1780,IM-1780T
Curved, No.7

Baby Jones Towel Clamp



IM-1781
Cross-Action, 66mm

Backhaus Towel Clamp



IM-1782
9cm

IM-1783
13cm

Serrefine Small



IM-1784
Straight, 4cm

IM-1785
Curved, 4cm

Hartmann Mosquito Forceps

95mm



IM-1786
Straight

IM-1787
Curved

Halstead

125mm




IM-1788
Straight

IM-1789
Curved

Jameson Muscle Forceps

With Slide Lock
6Teeth



IM-1790
Right

IM-1791
Left

Ptosis Lid Clamp



IM-1792
6-Pin, Concave

Berke Ptosis Forceps

With Slide Lock



IM-1793
20mm Jaws

IM-1794
27mm Jaws

Lambert Chalazion Forceps



IM-1795
Round, 8mm Dia.

Hunt Chalazion Forceps

Delicate 95mm



IM-1796
Round, 12mm Dia.

Desmarres Chalazion Forceps



IM-1797
Small, 20mm

IM-1798
Medium, 24mm

IM-1799
Large, 32mm

Snellen Entropion Forceps



IM-1800
Right

IM-1801
Left

Watzke Sleeve Spreading Forceps



IM-1802
Serrated 10mm tips, Angled 30°


Lens Holding Forceps



IM-1 803T

23Gauge, Delicate End, 90° Jaws, Squeeze Handle Mechanism


ICL Lens Cartridge Loading Forceps



IM-1 804T

20 Gauge, Straight

Micro Lens Cutter



IM-1 805T

For cutting soft acrylic and silicone IOLs
4mm long Stainlesssteel notched blades.

Lens Cutting & Holding Forceps



IM-1806T

To cut and remove all types of foldable intraocular Lens
Flat Titanium Handle Double Action, StainlessSteel Tips


Lens Cutting Forceps



IM-1 807T

Blade Length 7.7mm
Flat Titanium Double Action,Stainless Steel Tips


Livernois Pick-up & Holding Forceps



IM-1808,IM-1808T
for Soft IOL in One Motion

Livernois Pick-up & Holding Forceps

IOL Holding Forceps
Highly Polished Conical Jaws



IM-1809,IM-1809T
Straight

IM-1810,IM-1810T
Angled

Faulkner Lens Holding Forceps



IM-1811,IM-1811T
Smooth jaws, Rounded

Kelman-McPherson Lens Holding Forceps



IM-1812,IM-1812T
Cross-Action

Ernest-McDonald IOL Inserting Forceps



IM-1813,IM-1813T
5mm Jaws, Cross-Action

IM-1814,IM-1814T
6mm Jaws, Cross-Action

Spaleck Soft IOL Inserting Forceps

For Foldable Acrylic and Thin
Silicone Lenses



IM -1815, IM-1815T
Extra Delicate

Bipolar Forceps Jewelers Type



IM -1816
Pointed tips, Straight

Bipolar Forceps McPherson Type



IM-1817
Blunt tips, Angled

Bipolar Coaptation Forceps



IM -1818
Pointed Tips

Sinsky Lens Holding Forceps



IM-1820,IM-1820T
Delicate T-Shaped

Jaffe-Blaydes Lens Holding Forceps



IM-1821,IM-1821T
Very Delicate, Angled 30°, Narrow Smooth Jaws

Clayman Lens Holding Forceps



IM-1822,IM-1822T
Delicate Angled Jaws

Kratz Lens Holding Forceps



IM-1823,IM-1823T
Very Delicate , Curved, Narrow Smooth Jaws


Shepard Lens Holding Forceps



IM-1824,IM-1824T
With Serrated Lower Jaw

Castroviejo Needle holder

Delicate Jaws, 140mm Long



IM-1831
Straight, With Lock

IM-1832(25678-14)
Straight, Without Lock

IM-1833
Curved, With Lock

IM-1834(25679-14)
Curved, Without Lock

Castroviejo Needle holder


Standard Jaws, 140mm Long



IM-1835(25670-14)
Straight

IM-1836(25671-14)
Curved

SilMocks Needle holder




IM-1837(25612-15)
145mm 5.75"

Kalt Needle holder



IM-1838(25628-14)
Delicates mooth jaws, Straight With Lock

Birks Needle holder



IM-1839,IM-1839T
Extra Delicate, Without Lock

Drews Needle holder

Extra Delicate Jaws



IM-1840,IM-1840T Straight, With Lock	IM-1841,IM-1841T Straight, Without Lock	IM-1842,IM-1842T Curved, With Lock	IM-1843,IM-1843T Curved, Without Lock
--	---	--	---

Barraquer Needle holder

Short Model, Very Delicate Jaws 10cm



IM-1844,IM-1844T (25656-10) Straight, With Lock	IM-1845,IM-1845T (25650-10) Straight, Without Lock	IM-1846,IM-1846T (25657-10) Curved, With Lock	IM-1847,IM-1847T (25651-10) Curved, Without Lock
--	---	--	---

Barraquer Needle holder


Very Delicate Curved Jaws 13cm



IM-1848,IM-1848T With Lock	IM-1849,IM-1849T Without Lock
--------------------------------------	---

Troutman Needle holder

Delicate, Strongly Curved Jaws



IM-1850,IM-1850T With lock	IM-1851,IM-1851T Without Lock
--------------------------------------	---

McPherson Needle holder

Extra Delicate Jaws



IM-1852,IM-1852T
Straight, With Lock

IM-1853,IM-1853T
Straight, Without Lock

IM-1854,IM-1854T
Curved, With Lock

IM-1855,IM-1855T
Curved, Without Lock

Extra Fine Needle Holder




IM-1856T
Straight, Length of Jaws 6mm
Flat Titanium Handle Double Action
Stainless Steel Tips

IM-1857T
Curved, Length of Jaws 6mm
Flat Titanium Handle Double Action
Stainless Steel Tips

Barraquer Needle holder

Standard Curved Jaws 13cm



IM-1858,IM-1858T(25731-13)
With Lock

IM-1859,IM-1859T(25661-13)
Without Lock

Anis Needle holder

Extra Delicate Jaws



IM-1860,IM-1860T
Straight, With Lock

IM-1861,IM-1861T
Straight, Without Lock

IM-1862,IM-1862T
Curved, With Lock

IM-1863,IM-1863T
Curved, Without Lock

Corneal Trephine Blades &Handles



IM-1896
UniversalTrephineHandle

IM-1881 6.0mm	IM-1882 6.5mm	IM-1883 6.75mm	IM-1884 7.0mm	IM-1885 7.25mm
IM-1886 7.5mm	IM-1887 7.75mm	IM-1888 8.0mm	IM-1889 8.25mm	IM-1890 8.5mm
IM-1891 8.75mm	IM-1892 9.0mm	IM-1893 9.5mm	IM-1894 14.0mm	IM-1895 16.0mm

Schiotz Tonometer




IM-1900(44406-00)
autoclavable,
inclined scale

33

Sterilizing Tray

With Silicon Mat
(25cmx15cmx1.9cm)




IM-1901

Autoclaving, Steam & Ethylene Oxide Sterilizing
Silicone Mat, Protection for Instruments
Sterilization and Storage

Sterilizing Tray


With Two Silicon Inserts
(17cmx6.5cmx2cm)



IM-1902

Autoclaving, Steam & Ethylene Oxide Sterilizing
Pyramid Design Base for Better Steam Penetration
Sterilization and Storage

Extra Capsular Cataract Extraction Set/331



IM-1911

- | | |
|---|---|
| 1. IM-1602 Bishop-Harmon Anterior Chamber Irrigator Complete | 19. IM-1749 Bishop-Harmon Tissue Forceps StandIM |
| 2. IM-1781 Baby Jones Towel Clamp Cross-Action, 66mm | 20. IM-1684 McPherson Tying Forceps 5mm Long Smooth Tips Straight |
| 3. IM-1786 Hartmann Mosquito Forceps Straight | 21. IM-1685 McPherson Tying Forceps 5mm Long Smooth Tips Angled |
| 4. IM-1784 Serrefine Small Straight 4cm | 22. IM-1686 Kelman-McPherson Forceps, Delicate, 7.5mm Long Smooth Tips |
| 5. IM-1013 Lieberman Style Wire Speculum 15mm, K-wire, StandIM | 23. IM-1707 Ultrata Capsulorhexis Forceps Very Delicate Grasping Tips |
| 6. IM-1002 Barraquer Heavy Wire Speculum 14mm | 24. IM-1284 Sinskey II IOL Manipulating Hook Angled, 0.2mm, Blunt Tip |
| 7. IM-1166 Scalpel Handles 10cm | 25. IM-1290 Kuglen Iris Hook & Lens Manipulator Angled, Push-Pull Model |
| 8. IM-1424T Step Diamond Knife Tri-facet | 26. IM-1226 Barraquer Iris Spatula Tip Width 0.50mm Round, Smooth |
| 9. IM-1502 Castroviejo Corneoscleral Scissors Small blades, Right | 27. IM-1849 Barraquer Needle Holder without Lock 13cm |
| 10. IM-1503 Castroviejo Corneoscleral Scissors Small blades, Left | 28. IM-1835 Castroviejo Needle Holder Straight |
| 11. IM-1489 Westcott Tenotomy Scissors Blunt tips, Medium blades, Curved | 29. IM-1244 Lewis & SNELLEN Lens Loop Large |
| 12. IM-1488 Jaffe Stitch Scissors, pointed tips, Medium blades, Curved | 30. IM-1595 Anis Cortex Aspirating Cannula 0.4mm side port, 23 Gauge, Angled |
| 13. IM-1511 Vannas Capsulotomy Scissors Curved, Blade Length: 7mm | 31. IM-1562 Rycroft Air Injection Cannula 27 Gauge, Angle to Tip: 5mm |
| 14. IM-1520 Capsulotomy Scissors Stern-Gills, Angled forwalIM, 10mm | 32. IM-1573 McIntyre Anterior Chamber Cannula 26 Gauge, Angled |
| 15. IM-1097 Castroviejo Caliper Straight, 20mm Spread | 33. IM-1579 Jensen Posterior Capsule Polisher Sand Blasted Tip, 23 Gauge |
| 16. IM-1763 Bonaccolto Conjunctiva Forceps 1.2mm wide, Serrated Jaws | 34. IM-1590 Irrigating-Aspirating Cannula Simcoe |
| 17. IM-1663 Harms-Colibri Forceps Very Delicate, 1x2 Teeth, 0.12mm | 35. IM-1627 Brierley Chamber Maintainer 20 Gauge, Tapered Oval Tip |
| 18. IM-1668 Bonn Forceps Short Handle, 1x2 Teeth, 0.12mm | 36. IM-1901 Sterilizing Tray with Silicon Mat |

Extra Capsular Cataract Extraction using Phacoemulsification Set/332



IM-1912


- | | |
|---|---|
| 1. IM-1602 Bishop-Harmon Anterior Chamber Irrigator Complete | 20. IM-1690 Tennant Tying Forceps Extra Delicate, Smooth Tips Curved |
| 2. IM-1013 Lieberman Style Wire Speculum 15mm, K-wire, StandalIM | 21. IM-1707 Utrata Capsulorhexis Forceps Very Delicate Grasping Tips |
| 3. IM-1007 Wire Speculum Temporal Spring to wall Nasal Side | 22. IM-1686 Kelman-McPherson Forceps ,Delicate, 7.5mm Long Smooth Tips |
| 4. IM-1411T Diamond Knives Tri-facet FH-MDS 5, 1x1, 0x0, 20mm / 45° / 30° | 23. IM-1813 Ernest-McDonald IOL Inserting Forceps 5mm Jaws, Cross-Action |
| 5. IM-1415T Diamond Knives Clear Cornea 90° FH-CC 6, 0x3,0x0,20mm /90° abg | 24. IM-1808 Livernois Pick-up & Holding Forceps for Soft IOL in One Motion |
| 6. IM-1166 Scalpel Handles 10cm | 25. IM-1848 Barraquer Needle Holder with Lock 13cm |
| 7. IM-1786 Hartmann Mosquito Forceps Straight | 26. IM-1705 Dodick Nucleus Cracker Cross Action |
| 8. IM-1781 Baby Jones Towel Clamp Cross-Action, 66mm | 27. IM-1590 Irrigating-Aspirating Cannula Simcoe |
| 9. IM-1349 Rosen Phaco Splitter Wedge Shaped | 28. IM-1568 Binkhorst Irrigation Cannula 22 Gauge Right |
| 10. IM-1290 Kugel Iris Hook & Lens Manipulator Angled, Push-Pull Model | 29. IM-1567 Binkhorst Irrigation Cannula 22 Gauge Left |
| 11. IM-1297 Bechert Nucleus Rotator Angled 9mm, Blunt tip | 30. IM-1608 Nucleus Hydrodissection Cannula 25 Gauge, Flat Angled Tip 8mm |
| 12. IM-1300 Lester IOL Manipulator 0.25mm dia. tip, Angled | 31. IM-1562 Rycroft Air Injection Cannula 27 Gauge, Angle to Tip: 5mm |
| 13. IM-1489 Westcott Tenotomy Scissors Blunt tips, Medium blades, Curved | 32. IM-1579 Jensen Posterior Capsule Polisher Sand Blasted Tip, 23 Gauge |
| 14. IM-1097 Castroviejo Caliper Straight, 20mm Spread | 33. IM-1433 Fine-Thornton Swivel Fixation Ring 13mm Dia. Blunt Teeth |
| 15. IM-1520 Capsulotomy Scissors Stern-Gills, Angled for wallIM, 10mm | 34. IM-1817 Bipolar Forceps Jewelers Type Pointed tips, Straight |
| 16. IM-1662 Colibri Forceps Very Delicate, 1x2 Teeth, 0.12mm | 35. IM-1819 Bipolar ColIM Autoclave able |
| 17. IM-1674 Castroviejo Suturing Forceps 0.12mm | 36. IM-1511 Vannas Capsulotomy Scissors Curved, Blade Length: 7mm |
| 18. IM-1722 Fechtner Conjunctiva Forceps Ring Shaped Tip | 37. IM-1901 Sterilizing Tray with Silicon Mat |
| 19. IM-1689 Tennant Tying Forceps Extra Delicate, Smooth Tips Straight | |



IM-1913

1. **IM-1019** Lieberman Style Wire Speculum 15mm, V-shaped
2. **IM-1465** Maloney Intra-operative Keratometer Measurement of Astigmatism
3. **IM-1705** Dodick Nucleus Cracker Cross Action
4. **IM-1284** Sinskey II IOL Manipulating Hook Angled, 0.2mm Dia. Blunt Tip
5. **IM-1651T** Buratto Bi-Manual Aspiration Hand piece 21 Gauge, Port Size 0.35mm
6. **IM-1644T** Buratto Bi-Manual Irrigation Hand piece 21 Gauge, Curved
7. **IM-1415T** Diamond Knives, Clear Cornea 90°, 2.8mm
8. **IM-1411T** Diamond Knives Tri-facet FH-MDS 5, 1x1, 0x0, 20mm / 45° / 30°
9. **IM-1561** Rycroft Air Injection Cannula 30 Gauge, Angle to Tip: 5mm
10. **IM-1290** Kuglen Iris Hook & Lens Manipulator Angled, Push-Pull Model
11. **IM-1707** Utrata Capsulorhexis Forceps Triangular Grasping Tips
12. **IM-1097** Castroviejo Caliper Straight, 20mm Spread
13. **IM-1489** Westcott Tenotomy Scissors Blunt tips, Medium blades, Curved
14. **IM-1488** Jaffe Stitch Scissors Very sharp pointed tips, Medium blades, Curved
15. **IM-1609** Nucleus Hydrodissection Cannula 25 Gauge, Flat Angled Tip 8mm
16. **IM-1608** Nucleus Hydrodissection Cannula 25 Gauge, Flat Angled Tip 8mm
17. **IM-1721** Paddle Prechopper Forceps, Angled
18. **IM-1718T** Akahoshi Micro Prechopper Forceps Tips Height 1.3mm
19. **IM-1901** Sterilizing Tray with Silicon Mat


Micro Incision cataract Surgery (MICS) Set/334



IM-1914

1. **IM-1416T** Diamond Knives Tapered Spear MICS FH-CC / 90° Angled
2. **IM-1299** Lester IOL Manipulator, Straight
3. **IM-1538T** Alio MICS Scissors 23 Gauge, Curved Shaft
4. **IM-1713T** Micro Incision Capsulorhexis Forceps, 23 Gauge
5. **IM-1700** Giannetti MICS Capsulorhexis Forceps Ultra Delicate Shanks Tips
6. **IM-1705** Dodick Nucleus Cracker Cross Action
7. **IM-1646T** Buratto Bi-Manual irrigating Hand piece 23 Gauge, Curved
8. **IM-1650T** Buratto Bi-Manual Aspiration Hand piece 23 Gauge, Curved
9. **IM-1901** Sterilizing Tray with Silicon Mat

Toric OL Set/335




IM-1915

1. **IM-1438** Mendez Degree Gauge
2. **IM-1439** Bores Axis Marker
3. **IM-1446** Toric Marker Half-ring for mark of meridian 12 mm
4. **IM-1580** Jensen Posterior Capsule Polisher, 25 Gauge
5. **IM-1902** Sterilizing Tray with Two Silicon Inserts

38

Glaucoma Set/336




IM-1916

- | | |
|---|--|
| 1. IM-1602 Bishop-Harmon Anterior Chamber Irrigator | 12. IM-1662 Colibri Forceps, 1x2 Teeth, 0.12mm |
| 2. IM-1781 Baby Jones Towel Clamp | 13. IM-1749 Bishop-Harmon Tissue Forceps, StandalIM |
| 3. IM-1786 Hartmann Mosquito Forceps, Straight | 14. IM-1097 Castroviejo Caliper, Straight |
| 4. IM-1026 Schott Eye Speculum 16mm wire blades | 15. IM-1422T Diamond Knife with Glaucoma Surgery Step Knife |
| 5. IM-1002 Barraquer Heavy Wire Speculum | 16. IM-1166 Scalpel Handles, 10cm |
| 6. IM-1471 Eye Scissors Straight, StandalIM | 17. IM-1201 Kelly Punch |
| 7. IM-1489 Westcott Tenotomy Scissors Blunt tips, Curved | 18. IM-1849 Barraquer Needle Holder, Curved without Lock |
| 8. IM-1483 Barraquer Iris Scissors, Blunt tips | 19. IM-1562 Rycroft Air Injection Cannula, 27 Gauge |
| 9. IM-1516 Capsulotomy Scissors, Vannas Scissors, Curved | 20. IM-1135 Harms Trabeculotomy Probe with Suture Holes Left |
| 10. IM-1488 Jaffe Stitch Scissors, Curved | 21. IM-1136 Harms Trabeculotomy Probe with Suture Holes Right |
| 11. IM-1722 Fechtner Conjunctiva Forceps Ring Shaped Tip | 22. IM-1901 Sterilizing Tray with Silicon Mat |

39


Corneal Transplant Set/337



IM-1917

- | | | | |
|--------------------|--|--------------------|---|
| 1. IM-1602 | Bishop-Harmon Anterior Chamber Irrigator | 14. IM-1890 | Corneal Trephine Blades & Handles 8.5mm |
| 2. IM-1781 | Baby Jones Towel Clamp | 15. IM-1896 | Corneal Trephine Blades & Handles Universal Trephine Handle |
| 3. IM-1061 | Flieringa Scleral Fixation Ring, Set of eight | 16. IM-1667 | Polack Double Corneal Suturing Forceps |
| 4. IM-1002 | Barraquer Heavy Wire Speculum | 17. IM-1674 | Castroviejo Suturing Forceps, 0.12mm |
| 5. IM-1033 | Lieberman Eye Speculum Solid Blades, Adjustable | 18. IM-1675 | Castroviejo Suturing Forceps, 0.3mm |
| 6. IM-1786 | Hartmann Mosquito Forceps, Straight | 19. IM-1663 | Harms-Colibri Forceps |
| 7. IM-1784 | Serrefine Small, Straight | 20. IM-1689 | Tennant Tying Forceps, Straight |
| 8. IM-1498 | Troutman-Castroviejo Scissors with stop, Small blades, Right | 21. IM-1690 | Tennant Tying Forceps, Curved |
| 9. IM-1499 | Troutman-Castroviejo Scissors with stop, Small blades, Left | 22. IM-1848 | Barraquer Needle Holder with Lock |
| 10. IM-1489 | Westcott Tenotomy Scissors, Blunt tips, Curved | 23. IM-1833 | Castroviejo Needle Holder Curved, with Lock |
| 11. IM-1488 | Jaffe Stitch Scissors, Curved | 24. IM-1561 | Rycroft Air Injection Cannula, 30 Gauge |
| 12. IM-1886 | Corneal Trephine Blades & Handles 7.5mm | 25. IM-1573 | McIntyre Anterior Chamber Cannula, 26 Gauge |
| 13. IM-1888 | Corneal Trephine Blades & Handles 8.0mm | 26. IM-1901 | Sterilizing Tray with Silicon Mat |


Deep Anterior Lamellar Keratoplasty (DALK) Set/338



IM-1918

1. **IM-1191** DALK Lamellar Dissector Semi-sharp dissector
2. **IM-1392** DALK Corneal Dissector Triangular Shaft
3. **IM-1509** DALK Micro Scissors Left + Right Set
4. **IM-1560** DALK Air Injection Cannula
5. **IM-1901** Sterilizing Tray with Silicon Mat

Lasik Set/339



IM-1919

1. **IM-1013** Lieberman Wire Speculum, K-wire
2. **IM-1405** Rowen Combination LASIK Spatula
3. **IM-1399** Hersh IntraLase Flap Spatula
4. **IM-1401** Fukasaku LASIK Spatula
5. **IM-1403** Katzen LASIK Flap "Unzipper"
6. **IM-1409** Thorlakson Z-LASIK Flap Lifter
7. **IM-1435** Gimbel Stabilizing Ring
8. **IM-1727** Buratto LASIK Flap Forceps Disc Shaped, Serrated Jaws
9. **IM-1614** Buratto LASIK Irrigation Cannula 25 Gauge, 3 Ports
10. **IM-1039** Lieberman Aspirating Speculum 15mm
11. **IM-1901** Sterilizing Tray with Silicon Mat

42

Lid Set/340



IM-1920

1. **IM-1602** Bishop-Harmon Anterior Chamber Irrigator
2. **IM-1781** Baby Jones Towel Clamp
3. **IM-1789** Halstead Mosquito Forceps, Curved
4. **IM-1788** Halstead Mosquito Forceps, Straight
5. **IM-1474** Stevens Scissors, Straight
6. **IM-1471** Eye Scissors, Straight
7. **IM-1166** Scalpel Handles, 10cm
8. **IM-1489** Westcott Tenotomy Scissors, Curved
9. **IM-1488** Jaffe Stitch Scissors, Curved
10. **IM-1794** Berke Ptosis Forceps 27mm Jaws
11. **IM-1792** Ptosis Lid Clamp 6 Pins
12. **IM-1849** Barraquer Needle Holder without Lock
13. **IM-1066** Jaeger Lid Plate
14. **IM-1834** Castroviejo Needle Holder Curved, without Lock

15. **IM-1097** Castroviejo Caliper, Straight
16. **IM-1089** Desmarres Lid Retractor, Size-1
17. **IM-1090** Desmarres Lid Retractor, Size-2
18. **IM-1091** Desmarres Lid Retractor, Size-3
19. **IM-1250** Graefe Muscle Hook, Size-2
20. **IM-1255** Fixation Hook Double, Sharp, Small
21. **IM-1763** Bonaccolto Conjunctiva Forceps, Serrated Jaws
22. **IM-1749** Bishop-Harmon Tissue Forceps, StandalIM
23. **IM-1675** Castroviejo Suturing Forceps, 0.3mm
24. **IM-1817** Bipolar Forceps Blunt tips, Angled
25. **IM-1819** Bipolar ColM Autoclaveable
26. **IM-1096** Helveston Tissue Retractor, 11mm wide
27. **IM-1087** Knapp Lacrimal Retractor, Sharp
28. **IM-1901** Sterilizing Tray with Silicon Mat

43

Oculoplastic Set/341



IM-1921

- | | | | |
|--------------------|---|--------------------|---|
| 1. IM-1032 | Lancaster Eye Speculum 16mm | 14. IM-1792 | Ptosis Lid Clamp, 6 Pins |
| 2. IM-1092 | Schepens Orbital Retractor | 15. IM-1797 | Desmarres Chalazion Forceps, Medium |
| 3. IM-1090 | Desmarres Lid Retractor, Size-2 | 16. IM-1786 | Hartmann Mosquito Forceps, Straight |
| 4. IM-1066 | Jaeger Lid Plate | 17. IM-1835 | Castroviejo Needle Holder, Straight |
| 5. IM-1101 | Wright Fascia Needle | 18. IM-1838 | Kalt Needle Holder Delicate, Straight with Lock |
| 6. IM-1086 | Knapp Lacrimal Retractor 4-Prong, Blunt | 19. IM-1181 | Freer Periosteal Elevator |
| 7. IM-1269 | Meyerhoefer Chalazion Curette, Size-3 | 20. IM-1184 | West Lacrimal Chisel, Straight |
| 8. IM-1273 | Bunge Evisceration Spoon Large | 21. IM-1131 | Wilder Lacrimal Dilator Size-2 |
| 9. IM-1256 | Fixation Hook Double Sharp Large | 22. IM-1124 | Bowman Malleable Lacrimal Probe Size 3 & 4 |
| 10. IM-1250 | Graefe Muscle Hook, Size-2 | 23. IM-1134 | Pigtail Probe with Suture Holes |
| 11. IM-1474 | Stevens Scissor, Straight | 24. IM-1163 | Scalpel Handles Suitable for blades |
| 12. IM-1478 | Enucleation Scissor, Medium Curved | 25. IM-1901 | Sterilizing Tray with Silicon Mat |
| 13. IM-1756 | Tissue Forceps, Straight | | |

Chalazian Set/342




IM-1922

1. **IM-1066** Jaeger Lid Plate
2. **IM-1267** Meyerhoefer Chalazion Curette, Size-1
3. **IM-1268** Meyerhoefer Chalazion Curette, Size-2
4. **IM-1269** Meyerhoefer Chalazion Curette, Size-3
5. **IM-1270** Meyerhoefer Chalazion Curette, Size-4
6. **IM-1472** Eye Scissors Curved, StandalM
7. **IM-1753** Manhattan Fixation Forceps
8. **IM-1796** Hunt Chalazion Forceps
9. **IM-1798** Desmarres Chalazion Forceps, Medium
10. **IM-1834** Castroviejo Needle Holder Curved, without Lock
11. **IM-1166** Scalpel Handles, 10cm
12. **IM-1901** Sterilizing Tray with Silicon Mat

45

Strabismus Set/343




IM-1923

1. **IM-1033** Lieberman Eye Speculum Solid Blades, Adjustable
2. **IM-1787** Hartmann Mosquito Forceps, Curved
3. **IM-1784** Serrefine Small, Straight
4. **IM-1094** Helveston Tissue Retractor, 7mm wide
5. **IM-1095** Helveston Tissue Retractor, 9mm wide
6. **IM-1096** Helveston Tissue Retractor, 11mm wide
7. **IM-1249** Graefe Muscle Hook Size-1
8. **IM-1250** Graefe Muscle Hook Size-2
9. **IM-1251** Graefe Muscle Hook Size-3
10. **IM-1097** Castroviejo Caliper, Straight
11. **IM-1489** Westcott Tenotomy Scissors, Curved

12. **IM-1488**affe Stitch Scissors, Curved
13. **IM-1668** Bonn Forceps Short Handle, 1x2 Teeth, 0.12mm
14. **IM-1725** Moody Fixation Forceps, 1x2 Teeth with Lock Left
15. **IM-1726** Moody Fixation Forceps, 1x2 Teeth with Lock Right
16. **IM-1760** Graefe Fixation Forceps
17. **IM-1691** Harms Tying Forceps, Straight
18. **IM-1835** Castroviejo Needle Holder, Straight
19. **IM-1848** Barraquer Needle Holder with Lock
20. **IM-1602** Bishop-Harmon Anterior Chamber Irrigator
21. **IM-1166** Scalpel Handles, 10cm
22. **IM-1901** Sterilizing Tray with Silicon Mat

46


Lacrimal Set/344



IM-1924

- | | | | |
|--------------------|---|--------------------|--|
| 1. IM-1602 | Bishop-Harmon Anterior Chamber Irrigator | 16. IM-1185 | Feer Lacrimal Chisel, Curved |
| 2. IM-1782 | Backhaus Towel Clamp, 9cm | 17. IM-1134 | Pigtail Probe with Suture Holes |
| 3. IM-1788 | Halstead Mosquito Forceps, Straight | 18. IM-1137 | Crawford Retrieving Hook |
| 4. IM-1789 | Halstead Mosquito Forceps, Curved | 19. IM-1127 | Bowman Malleable Lacrimal Probe Set of 6 Sizes |
| 5. IM-1204 | Kerrison Rongeur Size-1, 3.0mm Wide | 20. IM-1754 | Dressing Forceps, Straight |
| 6. IM-1141 | Nasal Speculum, Adult | 21. IM-1752 | Sauer Fixation Forceps |
| 7. IM-1138 | Beyer Rongeur | 22. IM-1764 | Bonaccolto Conjunctiva Forceps, Serrated Jaws |
| 8. IM-1143 | Mallet Lead Filled | 23. IM-1675 | Castroviejo Suturing Forceps, 0.3mm |
| 9. IM-1166 | Scalpel Handles, 10cm | 24. IM-1086 | Knapp Lacrimal Retractor, Blunt |
| 10. IM-1081 | Stevenson Lacrimal Sac Retractor | 25. IM-1475 | Stevens Scissor, Curved |
| 11. IM-1622 | Lacrimal Cannula, 23 Gauge, Malleable Tip, Straight | 26. IM-1471 | Eye Scissors, Straight |
| 12. IM-1623 | Lacrimal Cannula, 23 Gauge, Malleable Tip, Curved | 27. IM-1833 | Castroviejo Needle Holder, Curved with Lock |
| 13. IM-1133 | Castroviejo Lacrimal Dilator, Double-Ended | 28. IM-1488 | Alte Stitch Scissors, Curved |
| 14. IM-1181 | Freer Periosteal Elevator | 29. IM-1508 | Tibolt Punctal Scissors |
| 15. IM-1184 | West Lacrimal Chisel, Straight | 30. IM-1901 | Sterilizing Tray with Silicon Mat |



Vitreous RetinalSet/345



IM-1925

- | | |
|--|--|
| 1. IM-1025 Castroviejo Eye Speculum, Adjustable | 19. IM-1674 Castroviejo Suturing Forceps, 0.12mm |
| 2. IM-1784 Serrefine Small, Straight | 20. IM-1675 Castroviejo Suturing Forceps, 0.3mm |
| 3. IM-1788 Halstead Mosquito Forceps, Straight | 21. IM-1749 Bishop-Harmon Tissue Forceps, StandalIM |
| 4. IM-1471 Eye Scissors, Straight | 22. IM-1833 Castroviejo Needle Holder Curved, with Lock |
| 5. IM-1602 Bishop-Harmon Anterior Chamber Irrigator Complete | 23. IM-1831 Castroviejo Needle Holder Straight, with Lock |
| 6. IM-1474 Stevens Scissors, Straight | 24. IM-1684 McPherson Tying Forceps, Straight |
| 7. IM-1475 Stevens Scissors, Curved | 25. IM-1094 Helveston Tissue Retractor, 7mm wide |
| 8. IM-1489 Westcott Tenotomy Scissors Blunt tips, Medium blades, Curved | 26. IM-1096 Helveston Tissue Retractor, 11mm wide |
| 9. IM-1488 Jaffe Stitch Scissors Very sharp pointed tips, Medium blades, Curved | 27. IM-1585 Gills Aspirating Cannula, 10mm Tip |
| 10. IM-1166 Scalpel Handles, 10cm | 28. IM-1819 Bipolar CoIM Autoclave able |
| 11. IM-1092 Schepens Orbital Retractor | 29. IM-1817 Bipolar Forceps, Blunt tips, Angled |
| 12. IM-1250 Graefe Muscle Hook, Size-2 | 30. IM-1100 Braunstein Caliper 3.5mm & 4mm |
| 13. IM-1247 Gass Retinal Detachment Hook | 31. IM-1527T Vitreo-Retinal Scissors, Curved (R 12mm) ,23 Gauge |
| 14. IM-1065 Schocket Scleral Depressor | 32. IM-1731T Vitreo-Retinal Forceps, Serrated Crocodile Forceps, 23 Gauge |
| 15. IM-1097 Castroviejo Caliper, Straight | 33. IM-1741T Vitreo-Retinal Forceps, Straight and End Gripping Jaws, 23 Gauge |
| 16. IM-1802 Watzke Sleeve Spreading Forceps | 34. IM-1745T Vitreo-Retinal Forceps, ILM Forceps, 23 Gauge |
| 17. IM-1763 Bonaccolto Conjunctiva Forceps | 35. IM-1901 Sterilizing Tray with Silicon Mat |
| 18. IM-1762 Nugent Utility Forceps, Serrated Jaws | |

Enucleation Set/346




10

IM-1926

1. **IM-1032** Lancaster Eye Speculum
2. **IM-1273** Bunge Evisceration Spoon Large
3. **IM-1272** Bunge Evisceration Spoon Small
4. **IM-1274** Wells Enucleation Spoon
5. **IM-1250** Graefe Muscle Hook, Size-2
6. **IM-1474** Stevens Scissor, Straight
7. **IM-1479** Enucleation Scissor, Strong Curve
8. **IM-1756** Tissue Forceps, Straight
9. **IM-1789** Halstead Mosquito Forceps, Curved
10. **IM-1901** Sterilizing Tray with Silicon Mat

Pterygium Set/347



IM-1927

1. **IM-1166** Scalpel Handle, 10cm
2. **IM-1032** Lancaster Eye Speculum
3. **IM-1684** McPherson Tying Forceps, Straight
4. **IM-1674** Castroviejo Suturing Forceps, 0.12mm
5. **IM-1189** Gill Corneal Knife, Straight
6. **IM-1834** Castroviejo Needle Holder Curved, without Lock
7. **IM-1471** Eye Scissors Straight, StandalIM
8. **IM-1763** Bonaccolto Conjunctiva Forceps, Serrated Jaws
9. **IM-1760** Graefe Fixation Forceps
10. **IM-1901** Sterilizing Tray with Silicon Mat



PRECISION
HIGH QUALITY
OPTIMUM
PERFORMANCE
ECONOMICAL

Contact Info

Mobile: 0046-724538624, 0092-3002125079

E-mail: admin@2smedeq.com

Web Address: www.2smedeq.com

Head Office

H # 59, Street 9, Sector G,
DHA Phase 5, Islamabad, Pakistan.

Branch Office

Tåtsplatsen 2, LGH 1523,
21110 Västra hamnen, Malmö, Sweden

Customer Service

Mobile:03362319529



WWW.2SLEVERANTORER.COM